
# HOWDY

senior

IT

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|  |    |   |           |
|--|----|---|-----------|
| Indications for use  | 27 |   |           |
| Warnings   | 27 |   |           |
| Further warnings and cautions  | 28 |   |           |
| General guidelines   | 28 | HOWDY SENIOR Maintenance                                |           |
| Device and accessories description   | 29 | Maintenance of the electronic unit                      | 44        |
| A) Sensored garment-HOWDY SENIOR textile unit  | 29 | Washing the textile unit and maintenance of the sensors | 44        |
|  | 30 | Withdrawn from service                                  | 45        |
| B) HOWDY SENIOR electronic unit  | 31 | Final warnings and recommendations                      | 46        |
| Packaging and Accessories  | 31 | <b>WARRANTY</b>   | <b>46</b> |
| Applicable Norms   | 34 | <b>ANNEX 1</b>  | <b>51</b> |
| How to use HOWDY SENIOR  | 34 |   |           |
| App installation on Smartphone/tablet  | 34 |   |           |
| Connection of the electronic unit to the Smartphone/tablet (1st installation)                | 34 |   |           |
| Assembly of HOWDY SENIOR   | 36 |   |           |
| Sign in and Starting the monitoring session  | 38 |   |           |
| How to leave the App   | 38 |   |           |
| Monitoring session   | 39 |   |           |
| Events notification System   | 39 |   |           |
| What to do in case of notification   | 41 |   |           |
| Main features  | 41 |   |           |
| Application messages   | 41 |   |           |
| What to do if the electronics fail to start  | 42 |   |           |
| What to do if the led light in the electronic unit flashes red while switching on the device | 42 |   |           |
| What to do in case of failed or wrong monitoring   | 42 |   |           |
| What to do in case of communication failure  | 42 |   |           |
| Switching off the electronic unit manually   | 43 |   |           |
| How to recharge the Electronic Unit  | 43 |   |           |
| Emergency shutdown of the electronic unit  | 43 |   |           |
| Internal battery life of the electronic unit   | 44 |   |           |

## INTRODUCTION

HOWDY SENIOR is a continuous and non-invasive monitoring device, capable to detect physiological parameters, such as heart rate, breath frequency, posture and movement. HOWDY SENIOR is aimed at adults, 16 years and over. In combination with the textile unit, HOWDY SENIOR detects, processes and transmits physiological parameters to a dedicated application.

HOWDY SENIOR is a medical device, based on technology developed and patented internationally by ComfTech srl; it offers reliability and high innovativeness guaranteeing at the same time comfort for the user.

Figure 1.

*Graphic representation of the textile units (on the left) and of HOWDY SENIOR (electronic unit on the right).*

## INDICATIONS FOR USE


Continuous and non-invasive monitoring of the physiological parameters such as heart rate, breathing rate, posture and movement for hospitals or home environment in adult patients, over 16 years of age. Temperature and blood saturation can be optionally monitored from third-party devices through the Howdy App. As a HOWDY SENIOR accessory, the app only communicates with other sensors when a HOWDY SENIOR is already connected.

## WARNINGS

Before using the HOWDY SENIOR device, carefully read the instructions contained in this manual and keep it for future reference.

•Please read all instructions carefully before using the device.

27

**SENIOR PRINTED  
TEXTILE UNIT  
CMF9TXM17**

**SENIOR SENSORIZED  
TEXTILE UNIT  
CMF9TXM11**

**WOMEN SENSORIZED  
TEXTILE UNIT  
CMF9TXM15**

**HOWDY SENIOR  
electronic unit**

- You must fully comply with all the cautions/warnings given on the device and in the user manual.
- HOWDY SENIOR, and in particular the textile unit, cannot be used and worn by anyone who is or might be allergic to cotton and/or silver fibre.
- HOWDY SENIOR is mainly made of soft textile materials. The wires don't have any sharp edges to prevent the occurrence of any contact injury. However, we recommend that you handle the electronic unit with care.
- HOWDY SENIOR is classified as a "BF" rated power supply (in accordance with the International Medical electric safety standard -Applied Part-).

It is recommended not to put into direct contact the snaps of the electronic unit and the electrodes of the Textile Unit neither with each-other nor with the mass/ earth.

•HOWDY SENIOR should not be worn by individuals with pacemaker

### FURTHER WARNINGS AND CAUTIONS

**WARNING:** No modification of the product is permitted

The use of different parts, such as cables, other than those specified for the HOWDY SENIOR Device its Textile Unit, may result in increased emissions or decreased immunity and in this case we cannot guarantee the correct operation of HOWDY SENIOR.

HOWDY SENIOR is based on wireless data transmission. Other electric and electronic devices could interfere with the signal transmission. At present, there is no evidence of radio interferences with other devices within the limits set by the HOWDY SENIOR Device. In any case, it is advised to keep all other devices at a minimum distance of 30 cm from the electronic unit or from the receiver.

### GENERAL GUIDELINES

• Do not leave your device within the reach of house pets and/or children.

28

• In case of sensors deterioration or loss of contact, the device might not operate correctly.

Furthermore:

• Based on data transmission via Bluetooth, the range of the product is about **10 metres in open field**, however some physical obstacles (such as walls, metallic barriers, furniture, objects in the house or in other premises) that stand between the electronic unit and the data receiver could impede the correct data transmission and reduce the range; for this reason, it is important to check the position of the user (for example sleeping in a prone position could cover the electronic unit) which could interfere with the receiver while transmitting the data.

The electronic unit has been designed and developed to be IP22-rated for added security of the users

In case of accidental or voluntary opening of the textile casing incorporating the electronic device, the product no longer has warranty of validity and proper operation. Therefore, it declines all responsibility for the reckless use and outside the warranty of validity.

**Charging:** The electronic unit is designed and manufactured in such a way as to **prevent and avoid operation during charging**. In any case, the device is charged through a NFC-enabled battery charger Qi certified (not supplied with the device).

Important Information about the built-in rechargeable

batteries

Rechargeable batteries shall not be exposed to direct sunlight, flame or high temperatures due to a possible risk of explosion of the same when exposed to such agents.

### **Electromagnetic compatibility (EMF)**

This device complies with all the electromagnetic compatibility standards (EMC). If used appropriately following the instructions in this manual, the device allows safe use as confirmed by the currently available scientific results.

**Water and humidity** : do not use the device next to water sources, such as bath tubs and washbasins or in a humid environment.

Under no circumstances should the electronic unit be immersed in water. If this is the case, the correct restoration of the operation of the same electronic unit is not guaranteed.

### **Heat**

Keep the device away from heat sources, such as radiators, stoves or other heat generating devices.

## **DEVICE AND ACCESSORIES DESCRIPTION**

HOWDY SENIOR combined with the textile unit makes up a system, consisting of:

- a) **Textile unit**: garment with integrated textile sensors
- b) **HOWDY SENIOR Electronic unit**: data collection and transmission hardware

## **29**

### **A) Sensored garment –HOWDY SENIOR textile unit**

The sensed garment combines technology, comfort, practicality and simplicity of use. The textile unit has the following features:

- 1– made of 80% Polyamide and 20% Elastane, with Oekotex certified yarns;
- 2– simple and comfortable design;
- 3– no need for any special maintenance - washing resistant, if complying with the washing indications provided in this manual;
- 4– integrated non-allergenic textile sensors, for detection of monoderivation electrocardiograph, heart rate, breathing signal and breathing rate;
- 5– the connection of the sensors to the electronic unit is guaranteed by textile components, therefore no rigid element or electrical wire is inserted into the garment;
- 6 – the connection between the electronic unit and the textile unit is made through simple nickel-free snaps.


Figure 2 shows the textile sensors embedded into the garments. The textile unit CMF9TXM17 share the same characteristics (3 textile ECG electrodes and one strain gauge)

## KEY

1 Textile sensors


2 Connection clip with electronic unit


30

Figure 2.

*Textile units of HOWDY SENIOR*

## B) HOWDY SENIOR electronic unit

## KEY

1 On/off button

2 LED status indicator

3 Nickel-free snaps

Figure 3.

*HOWDY SENIOR electronic unit*

## PACKAGING AND ACCESSORIES

### HOWDY SENIOR System Supply equipment

- 1 textile unit
- 1 electronic unit
- 1 Use and Maintenance Manual

The packaging for transport consists of a cardboard box.

### HOWDY SENIOR Supply equipment

- 1 electronic unit
- 1 Use and Maintenance Manual

The packaging for transport consists of a cardboard box.

### HOWDY SENIOR Textile Unit Supply equipment

- 1 textile unit
- 1 Use and Maintenance Manual

The packaging for transport consists of a polybag.

#### Note:

To charge the electronic unit it is necessary to get an induction battery charger (NOT supplied) Qi standard compliant. The Qi standard identifying symbol must be displayed as shown below.

#### Note: Is Bluetooth technology dangerous to human health?

According to current research, the Specific Absorption Rate (SAR) is the basic factor to measure possible harm to our health caused by the exposure to high frequency electromagnetic fields.

The Absorption Rate describes the amount of RF energy (measured in Watt (W), 1 W= 1000 mW) absorbed by the human body (Kg).

**31** The maximum SAR allowed is 0,08 W/Kg for the whole body and 2,00 W/Kg for body parts, for example for the head.

The devices connected via Bluetooth or WLAN, as single components, remain far below this SAR limit. HOWDY SENIOR has been tested for electromagnetic compatibility and totally complies with the SAR limit. When the limit is respected, at the current state of science, there is no evidence that high frequency electromagnetic fields can be a risk to one's health.


## APPLICABLE NORMS

HOWDY SENIOR is classified as a class IIb medical device, in accordance with the requirements of 93/42/CEE Directive on Medical Devices amended by the 2007/47/CE Directive.

The HOWDY SENIOR meets all the requirements of the following regulations:


- IEC 60601-1
- IEC 60601-1-2
- IEC 60601-2-47
- EN ISO 15223-1
- EN 1041
- EN ISO 14971
- EN ISO 10993-1
- EN 60601-1-6
- EN 62304
- EN 62366
- EN 60601-1-11
- ETSI EN 301 489-17 V3.1.1
- EN 301 489-3 V2.1.1
- ETSI EN 301 489-1 V2.1.1
- EN 62479:2010
- EN 300 328 V2.1.1


**LIST OF USED SYMBOLS – HOWDY  
SENIOR ELECTRONIC UNIT**

32

**LIST OF USED SYMBOLS – TEXTILE UNIT  
(see label)**

|  |  |   |   |
|--|--|---|---|
|  | Manufacturer                           |  | Read the Instructions prior to use                        |
|  | Attention                              |  | CE Mark   |
|  | Product code                           |  | Do not dispose as household waste - it contains a battery |
|  | Lot Number                             |  | Type BF applied part – isolated                           |
|  | Do not use if the casing is not intact |  | Induction NFC Recharge                                    |
|  | Pressure limit: 70kPa – 106 kPa        |  | Charge  |
|  | Keep away from humidity                |  | Wireless Device complying to Qi standard                  |
|  | Humidity limit: 10% - 95%              |  | Bluetooth   |
|  | IP Protective Rating                   |  | Temperatura di funzionamento tra -10°C e 70°C             |
|  | Non-ionising Radiations (NIR)          |  | Contains Lithium-ion battery                              |
|  | CE Mark and Notified Body Number       |  | Textile unit size indication                              |

## HOW TO USE HOWDY SENIOR

HOWDY SENIOR works in combination with a smartphone/tablet that has sufficient battery to guarantee a Bluetooth transmission.

It is recommended to check this every time you want to start monitoring. In alternative, you can use HOWDY SENIOR with a smartphone/tablet connected to the electricity supply.

The battery charge level is displayed in the main screen of the dedicated app and you can also get a notification by a flashing light in the electronic unit, as explained in the following paragraphs.

HOWDY SENIOR can be used both by qualified healthcare professionals and by adults who do not operate in the healthcare sector. In the latter case, it is required a minimum digital competence, an adequate comprehension of the operational instructions and absence of severe vision or hearing problems. If the operator is also a patient, all the functions listed below can be used.

### APP Installation on Smartphone/tablet

HOWDY SENIOR works in combination with a dedicated app that displays the data monitored by HOWDY SENIOR.

The recommended operating system of the smartphone/tablet should be Android 8.0 or above. Follow the installation wizard emphasizing the required consents during the procedure

Before connecting the Electronic Unit to the Smartphone/Tablet it is necessary to activate its

### Connection of the electronic unit to the Smartphone/tablet (1st installation)

The electronic unit transmits the patient's heart rate and movement data to the smartphone/tablet via a Bluetooth Low Energy wireless connection. Therefore, first of all it is necessary to switch on the Bluetooth signal and the position/location of the smartphone/tablet. Check the instructions manual to find out how to switch on the Bluetooth signal and the position/location of the smartphone/tablet that you are using.

The first time you use it, you need to select the Electronic Unit you intend to connect to from the Options Menu of the App: touch on "Add Sensor" and then choose the Electronic Unit you wish to use (refer to the section "Sign in and Starting the monitoring session"). You don't need to repeat this operation the next times you use the app; the connection will be automatic.

### ASSEMBLY OF HOWDY SENIOR

**1) Switch on the electronic unit pressing the small button marked by the ignition switch symbol for at least two seconds.** Every time the electronic unit is switched on any initial check is done on the main hardware components. This phase is identify by a green LED followed by a series of red blinks. If the initial check is positive, it will start the power-on sequence: green LED-Red blinking LED- blue blinking LED, which indicates that the device is ready for the Bluetooth connection (figure 5).


Figure 5.

*Switching on HOWDY SENIOR electronic unit*


**2) Connect the electronic unit to the textile unit by pressing the snap buttons together.** Pay attention to the correct direction of assembly: the dash of the ignition switch symbol  must be positioned towards the face as shown.

a. We recommend to exert a moderate pression, on one side first and then on the other side of the unit, while putting the connection clips together.

**Note:** For a correct electronic connection both automatic buttons need to be well in contact with each-other (when you connect the snaps correctly you will hear a distinctive fastening clicking noise).

It is also recommended to keep the snaps clean, clear of dust and other agents, in order to have a good electric connection between the metallic terminals.

**3) How to wear or make someone wear a garment (textile unit).** When using the SENIOR SENSORIZED TEXTILE UNIT CMF9TXM11, put the top on after unfastening the side fastener and adjust it so that the sensors fit the chest. When using the WOMEN SENSORIZED TEXTILE UNIT CMF9TXM15 or the SENIOR SENSORIZED PRINTED BELT MF9TXM17, put the garment on after opening the Velcro fastener and adjust the fastener so that the sensors fit snugly against the chest.

The electronic unit will remain covered by the fabric band in the first case.

**4) Make sure that the textile sensors inside the garment are well in contact with the body;**  
o you can test it inserting a finger inside the garment and smoothing the sensors out, avoiding creases and turn-ups;

o If a sensor rolls up or loses contact with the body, you just need to put the sensor in contact with the skin again and wait for a few seconds;

**Note:** the sensors are more effective if the skin is clean; grease and dust in the contact area could interfere with the radio signal.

**Note:** the sensors are more effective if slightly moisten with water. Sweat could also be an adjuvant for the quality of the electrode-skin interface.

**5) Start monitoring (as per instructions in the next paragraph)**

**Note:** After removing the electronic unit, put it on the charging mat and check that it is charging correctly (LED light on the charging unit will switch ON).

**Note:** The textile unit can be comfortably worn even without the electronic unit (in this case the monitoring will not work) to be ready when necessary.

## SIGN IN AND STARTING THE MONITORING SESSION

Sign in via the “Start” button. For communication with an external portal, log in using the credentials provided by the provider of the web portal. In this way the user’s data are sent to the portal.

Figure 6.  
Log-in Screen

- 1) Check that the Bluetooth and the Position/Location of your smartphone/tablet are turned on.
- 2) Switch on the electronic unit by lightly pressing on the ON/OFF button; a red LED light will start flashing. When the device is ready a blue LED light will start flashing.
- 3) Activate the App by touching with your finger the corresponding icon. As soon as an external device connects to the Electronic Unit, the blue LED light

36

will turn off and the “Sensor connected” notification will be displayed on your smartphone/tablet. On your smartphone/tablet display you will see this start screen:


Figure 7.

### Available Functions

*The Homepage is identified upwards and shows the parameters measured by active sensors. The footer shows report functions and the mask for browse between Home, Options, Notifications.*

The first time you use it, you need to select the Electronic Unit you intend to connect to on the smartphone/tablet from the App Settings menu: choose

“Add Sensor” and then the Electronic Unit you wish to use (HOWDY SENIOR is identified by BERTA\_x.x.x). The app will ask to match the added device with the type of corresponding sensors (BERTA\_x.x.x shall be associated with HOWDY SENIOR).

As soon as an external device (smartphone/tablet) connects to the Electronic Unit, the blue LED light will turn off and the “Sensor connected” notification will be displayed on your smartphone/tablet.

On the start screen you can now view the following parameters: heart rate, breathing rate, posture .

The Howdy App allows communication with additional sensors:

- Howdy Breath/Howdy Breath-1 : Manufacturer: ComfTech srl

- Thermometer:Tempsitter Wireless Thermometer; European Authorized Representative : Shangai International Holding Corp. GmbH(Europe)

- Wireless pulse oximeter : OXY-10 Pulse oximeter; Manufacturer: GIMA

Follow the procedure described above to connect the last two sensors to the app.

The thermometer will be displayed as: BLT\_MODT

The pulse oximeter will be displayed as: PC\_60NW

These sensors can communicate with the app if the Howdy device is already connected.

On the top right of the display it is shown the remaining charge of the Electronic Unit. We recommend to start using the device with a minimum of 20% charge. Please refer to the section **“HOW TO RECHARGE THE ELECTRONIC UNIT”** to recharge the device.

4) Through the “Generate report” button it is possible to obtain a pdf file with summary data of the monitoring. Please refer to the paragraph “MONITORING

SESSION” for more information.

5) Through the “Real-time ECG” button or clicking the parameters “Heart rate” and “Breath rate” button it is possible to view the graphs.

WARNING: the ECG tracing viewed through the APP has only qualitative value and it should not be intended and/or used for diagnostic purposes.

6) You can leave the APP at any time, following the instructions in the paragraph “HOW TO LEAVE THE APP”.


Figura 8.  
Settings  
Adding and managing sensors.

## HOW TO LEAVE THE APP

To leave the App, you can use the following buttons of your smartphone/tablet as described below.

“Home” button: the acquisition screen is hidden but it keeps running.

“Back/Return” button: press twice to leave the app. We recommend to close the App only once the monitoring session has ended.

## MONITORING SESSION

The monitoring session starts automatically with the starting of the device. The app keeps the data of the last 48 hours. With the function “Generate report” it is possible to select the period of interest among the last 48 hours and visualize the monitored data (also in case of interruptions of the monitoring session).

38

While monitoring the physiological parameters measured by the sensors will be updated in real time. You can generate a report at any time. With the function generate report it is possible to enter the operator ID if necessary and select the period of interest among the last 48 hours. It is then possible to visualize the report confirming the popup which appears.

The monitoring report consists of three parts: a header, a generic table, with the aggregated physiological data for the whole monitoring period, and some periodic tables with the physiological data divided into time periods (according to the duration of the monitoring). With the “Report Archive” button you can access all the generated reports and visualize them, send them via e-mail, or delete them swiping to the right and pressing the Trash Can symbol.

**Note:** ComfTech srl does not take any responsibility about the use of such files.

WARNING: the tracings recorded through the APP have only not be intended and/or for diagnostic purposes.


Figure 9.  
Report generation screen

Confidential

Confitech srl declines any responsibility for using this report document outside the Howdy app

18 March, 2021

|                     |   |
|---------------------|---|
| Monitored patient   |   |
| Operator            |   |
| Software version    | 1.1.0   |
| Start of monitoring | 18/03/2021 00:00                                  |
| End of monitoring   | 18/03/2021 07:40                                  |
| Monitoring time     | 7 Hours 40 Minutes                                |
| Device              | HOWDY SENIOR (fw v. 2.0.58) - (D5:A4:7E:62:88:E0) |

Figura 10.

*Header of monitoring report*

## EVENTS NOTIFICATION SYSTEM

Notifications can be visualized pressing “Notifications” button.

The notified events will then be summarized in the final reportfile.

The events can be of a technical nature, identified with yellow (eg: low battery in the Electronic Unit, device disconnected) or of a physiological nature, identified by red, such as heart rate/breathing rate outside the configured thresholds..

Default thresholds are: Heart rate: 50-120 bpm

Respiratory frequency: 8-24 rpm

The “Modify Thresholds” button allows, by accessing with the admin password, to change the Thresholds of the physiological parameter and to select the parameters of whom you want to receive notifications.

## WHAT TO DO IN CASE OF NOTIFICATION

**a)** Go immediately to the user to check his/her health condition. If the user is also the operator, go straight to point **b)**

**b)** If the user has no obvious health problems, check the following possible conditions:

- Check if the electronic unit has detached from the textile unit.
- Check if the textile sensors are in the correct position and are in direct contact with the skin of the user.
- Check if the electronic unit is switched on.
- Check if the battery of the electronic unit has sufficient charge.
- Check if the Bluetooth connection and the data transmission is on.

**c)** If none of the above conditions are met, therefore you believe that the device is working correctly, we recommend to call the doctor.



Figure 11.  
Screen menus of Notification Settings

The main features of the device are the following:

- 1) Measuring the heart rate;
- 2) Measuring the breathing rate;
- 3) Data transmission of heart rate and breathing rate via Bluetooth to a dedicated smartphone/tablet .

## APPLICATION MESSAGES

| Message  | Possible Causes   | Actions  |
|--|---|--|
| Low battery  | The battery level of the electronic unit is below 20%   | It is recommended to stop all activities and put the electronic unit on the charging mat.  |
| Sensor Connected                                   | The connection between the electronic unit and the smartphone/tablet is active                | HOWDY SENIOR is ready for use  |
| Sensor Disconnected                                | There is NO connection between the electronic unit and the smartphone/ tablet                 | 1. Select once more the desired device from the App settings<br>2. If you get the same message again, restart the electronic unit and the App. |
| Heart rate above...<br>Heart rate below...         | The heart rate detected by HOWDY SENIOR is outside the normality range or values not received | Check the paragraph WHAT TO DO IN CASE OF NOTIFICATION on page 39  |
| Breathing rate above...<br>Breathing rate below... | The breathing rate detected by HOWDY SENIOR is outside the normality range                    | Check the paragraph WHAT TO DO IN CASE OF NOTIFICATION on page 39  |

## WHAT TO DO IF THE ELECTRONICS FAIL TO START

If the electronic unit fails to start, it is very likely that the battery is out of charge. We recommend to recharge it for at least 30 minutes and then try to start it again.

## WHAT TO DO IF THE LED LIGHT IN THE ELECTRONIC UNIT FLASHES RED WHILE SWITCHING ON THE DEVICE

Wait for about one minute and check that the LED flashing light turns to blue.

If that doesn't happen, switch off the device pressing for about two seconds on the power button. The LED light should switch off straight after. Repeat the procedure to switch on the device.

If the problem still persists, please contact the Manufacturer. You will find the contact details on the first page of this manual.

## WHAT TO DO IN CASE OF FAILED OR WRONG MONITORING

1. First of all, check the health condition of the user. If the user is also the operator, go straight to point 2.

2. If the user is ok and the situation remains the same, check that the electronic unit is correctly connected to the snaps of the textile unit;

3. In case this unexpected situation persists, check that the textile electrodes are correctly in contact with the skin of the user;

4. In case this unexpected situation persists, the textile unit might be damaged; we recommend to check the state of the textile sensors in order to decide whether to purchase a new textile unit.

Please be aware that during the first 3 minutes of monitoring, the system needs to adapt to the unique morphology of the user.

## WHAT TO DO IN CASE OF COMMUNICATION FAILURE

42

If the Electronic Unit does not connect to the device, it is possible that it is not correctly paired to the smartphone/tablet where the App is installed. We recommend to repeat the pairing procedure as described on page 34 of this user manual.

In case the Electronic Unit were correctly paired, but the App is not able to establish the proper connection, we recommend to switch off the Electronic Unit and turn it on again. Likewise, close and reactivate the App before trying again. Also, please check that the position/location on the smartphone/tablet is on.

**Note:** if, after following the steps listed above, there is still no communication between the electronic unit and the smartphone/tablet where the App is installed, please contact Customer Support – you will find the details on the first page of this manual.

## SWITCHING OFF THE ELECTRONIC UNIT MANUALLY

The electronic Unit can be switched off manually at any time.

To do so, detach the electronic unit from the textile unit and press the ON/OFF button for at least 3 seconds. This lapse of time is necessary to prevent the accidental shutdown of the electronic unit by the user. After shutting down, the LED light in the electronic unit will flash red and the same LED light will then switch off. In case this operation is done while the connection between App and electronic unit is active, you will receive a "Disconnected" notification message on your smartphone/tablet's display to alert you that the electronic unit has been disconnected

## HOW TO RECHARGE THE ELECTRONIC UNIT

To recharge the electronic unit, position it with its buttons facing down on the induction charging mat (NOT supplied with the device).

During charging, all the electronic parts (except those ones necessary for recharging the battery) are automatically switched off for the safety of the user. It is important to verify that the device is charging checking that the LED light on the charging mat is ON.

The electronic unit must be positioned in correspondence with the charging mat's antenna.

When you position the electronic unit on the charging mat, the LED light will flash red and green and it will flash only green when the charging is complete.

**FLASHING GREEN and RED:** the electronic unit is charging

**FLASHING GREEN slowly:** the electronic unit is charging and is over >95% full

**FLASHING GREEN:** the charging of the electronic unit is complete

## EMERGENCY SHUTDOWN OF THE ELECTRONIC UNIT

At any time, the electronic unit can be shut down as an emergency and this will guarantee a resetting of the unit, also when it has entered a state of error which cannot be fixed during normal operation. To switch it off, keep pressing the button for about 15 seconds.

## INTERNAL BATTERY LIFE OF THE ELECTRONIC UNIT

The electronic unit contains an internal LiPo battery pack with a capacity of 600 mAh rechargeable with NFC technology (battery life 30 hours used with continuous

data transmission through Bluetooth® connection). Normally, this type of batteries have a satisfactory charging capacity of 500 refills, equal to about 2 years of life.

It is recommended to review weekly the length of monitoring done: if, after full charging, the battery cannot guarantee at least 8 hours of continued data transmission via Bluetooth, it is advised to replace the electronic unit with a new one.

**DO NOT** open the electronic unit and **DO NOT** try to substitute the internal battery as the electronic unit would be irretrievably damaged!!!

Warning to resume charging after a long period of non-use of the device:

If the electronic unit hasn't been used for over a month, before using it again we recommend to execute at least 2 charge and discharge cycles of the electronic unit itself.

In case of non-compliance with this rule and the prior verification of the battery life is not liable for damage to property and people.

## CHECKING THE OPERATION STATUS

The electronic unit, if it's not connected to any smartphone/tablet, will send an alert message with the flashing blue text: "ready for connection".

**FLASHING BLUE:** the electronic unit is ready for connection with the smartphone/tablet.

Once connected, the electronic will no longer flash not to bother the user or anyone around.

To check the operation status of the electronic unit,

press the button briefly:

**FLASHING GREEN ONCE:** the electronic unit is transmitting and is charged (>20%)

**FLASHING GREEN TWICE:** the electronic unit is saving the recording in the memory and is charged (>20%)

**FLASHING RED ONCE:** the electronic unit is transmitting and its charge is low (<20%)

**FLASHING RED TWICE:** the electronic unit is saving the recording in the memory and its charge is low (<20%)

## MAINTENANCE OF HOWDY SENIOR

The maintenance of HOWDY SENIOR can be done both by healthcare professionals and by adults who do not operate in the healthcare sector. In the latter case, it is required a minimum digital competence, an adequate comprehension of the operational instructions and absence of severe vision or hearing problems. If the operator is also a patient, all the functions listed below can be used.

To make sure that the HOWDY SENIOR Device is operating correctly, we recommend to carry out the following operational checks at least once a week:

1. Visual inspection of the hermetic closure of the electronic unit: there should be no openings.
2. Visual inspection of the textile unit: the sensors should be smooth and with an even surface.

## 44

### Environmental conditions for use, transport and storage of the HOWDY SENIOR Device

- Operating temperature range: +5°C – +45 °C
  - Operating relative humidity: 10% - 95% without condensation
  - Atmospheric operating pressure: 700-1060 hPa
  - Temperature range for storage/depot: -10 – +70 °C
  - Relative humidity in storage/depot: 10% - 95% without condensation
  - Atmospheric pressure in storage/depot: 700-1060 hPa
- 
- Protection level of the electronic unit: IP22 (Protected against solid bodies of superior dimensions to 12mm, protected against falling water when tilted up to 15°).

### MAINTENANCE OF THE ELECTRONIC UNIT

- Please use a cloth damped with water to clean the electronic unit.
- Please take care when using the damped cloth and avoid water or liquid next to the snap buttons connecting the electronic unit to the textile unit.

The electronic unit cannot be submerged in water: always remember to detach it from the textile unit before washing it.

### WASHING THE TEXTILE UNIT AND MAINTENANCE OF THE SENSORS

- It is recommended to wash the textile unit after each use.
- The textile unit should normally be washed in the washing machine at 30° with delicate cycle.

For hygiene reasons, the product is to be considered personal. In any case, after about 6 months of use, we suggest to purchase new textile units.

### WITHDRAWN FROM SERVICE

The textile unit (except for the sensorised belt CMF9TXM17) is a normal textile product which can be disposed of according to the standard procedures set for this type of product.

However, the electronic unit of HOWDY SENIOR and the sensorised belt CMF9TXM17 need to be disposed of according to the RAEE Directive guidelines, therefore once discarded, they should be disposed of through a separate refuse collection and delivered to an “ecological island” set up by the responsible municipality or by the equipment distributor, even when purchasing a new one.

With the “RAEE Directive” Italy approved the European Directive 2012/19/EU about the use of particular hazardous substances in Electric and Electronic Equipment and about managing their disposal.

To protect the environment, and send these products (RAEE) to the adequate waste centre, it's required a differentiated collection of this type of equipments. At present, all citizens can send their own waste to ecological islands; since 18th June 2010, it is possible to return the waste directly to the seller free of charge, when purchasing a similar equipment. This procedure is regulated by the decree D.M.65 of 8th March 2010. The distributors (individuals who sell new products to end users) are obliged to withdraw the equipment to be discarded at the moment of purchase of a new, equivalent equipment. Illegal disposal of the above-mentioned products shall be subject to administrative

- Do not soak it..
- Avoid any abrasive products or with chemical solvents. It is advised to use a liquid detergent.
- The Manufacturer declines all liability for damage to the textile unit due to the use of any fabric conditioner, as this causes a rapid deterioration of the textile electrodes, totally impairing the sensing function.
- Do not bleach the textile unit and do not submerge it in irritating disinfectants (in case this is absolutely necessary, avoid contact of these substances with the textile sensors).
- The textile unit can be lightly spun.
- The textile unit cannot be put in a tumble-dryer: this causes a rapid deterioration of the textile electrodes, totally impairing the sensing function.
- Because of the particular weave, the textile unit does not need to be ironed; in any case, always avoid a direct contact of the iron with the textile sensors.

These instructions are reproduced on the garment label.

The textile sensors are the most “sensitive” part of the textile unit, therefore the most delicate:

- Always keep them inside the textile unit; when drying do not expose to direct sunlight (it is sufficient to dry the garment without overturning it).
- Do not put the iron over the sensors.
- Do not twist them when wringing by hand.

If you follow the maintenance instructions listed above, your textile unit will maintain its functions unaltered and will maintain its sensitive capabilities even after repeated washing (up to 30 washings).

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The symbol RAEE is shown in the table at page 8.

## FINAL WARNINGS AND RECOMMENDATIONS

Under no circumstances shall the electronic unit be tampered with, as this may alter the correct functioning of the device and in this case any product and performance warranty will be rendered invalid. Should the product need to be repaired, this can be done only by qualified personnel.

For any technical problem, please contact Customer Support - refer to the manufacturer's contact information provided at the end of this manual.

On the outer packaging of HOWDY SENIOR, you can find: the serial number unique to the product itself, the product certification, technical information and the manufacturer's details according to regulations (also listed at the end of this manual).

For any problem or further information please contact

**ComfTech S.r.l.**  
P.za Castello 9,  
20900 Monza (MB)  
Tel. +39 039 9008300  
[www.comftech.com](http://www.comftech.com)

HOWDY SENIOR is covered by a 24 months guarantee, from the date of purchasing, against any material or manufacturing defects.

The 24 months warranty, pursuant to Legislative Decree n. 24/2002, applies to products with a lack of conformity, provided that it is used correctly, in compliance with its intended use and with the provisions of the attached technical documentation. This warranty, in compliance with Legislative Decree no. 24/92, is reserved only to private consumers (natural person who purchases the goods for purposes not related to his/her professional activity, or makes the purchase without indicating in the order form a reference VAT number). In case of lack of conformity, ComfTech s.r.l. provides, at no cost to the Customer, the restoration of product conformity by repair/ replacement, until termination of the contract.

**Important:** Immediately fill out the purchase coupon below and attach the receipt to the coupon. Keep the documentation in a safe place.

If you have a complaint, please contact the distributor from which you purchased the system directly. Remember to keep the receipt; you may need it as proof of purchase.

The manufacturer will examine the device and, if it is found to be defective due to materials or incorrect processing, will replace the device. If the device covered by the warranty has been damaged due to misuse, has been modified or repaired by an unauthorised person or has been damaged, the warranty will not be valid. The warranty is personal and non-transferable.

**The operation of the textile unit is guaranteed only and exclusively with its electronic unit. We decline any responsibility in the event that the textile unit is used with devices other than the electronic unit.**

EXCEPT AS INDICATED ABOVE, THE PRODUCT IS SOLD "AS IS" WITHOUT ANY EXPRESS OR IMPLIED WARRANTIES FOR THE PRODUCT WHATSOEVER, INCLUDING, WITHOUT LIMITATION, ANY IMPLIED WARRANTIES REGARDING THE CONDITION, DESIGN, TECHNICAL SPECIFICATIONS, WORKMANSHIP, MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OF THE PRODUCT, OR ANY WARRANTIES THAT THE PRODUCT IS FREE FROM LATENT DEFECTS OR DEFICIENCIES, OR THAT THE PRODUCT IS FREE FROM INFRINGEMENT OF ANY PATENT, TRADEMARK, COPYRIGHT OR PROPRIETARY RIGHT OF ANY THIRD PARTY AND SUCH WARRANTIES ARE HEREBY EXPRESSLY DISCLAIMED.

THIS LIMITED WARRANTY APPLIES ONLY AND EXCLUSIVELY FOR ALLEGED PRODUCT FAILURE. IN NO EVENT SHALL WE BE LIABLE FOR ANY INDIRECT CONSEQUENTIAL, SPECIAL, ACCIDENTAL OR PUNITIVE DAMAGES, WHETHER IN AN ACTION IN CONTRACT, TORT, STRICT LIABILITY, OR IMPOSED BY STATUTE, OR OTHERWISE. OUR LIABILITY SHALL IN NO EVENT EXCEED THE PURCHASE PRICE OF THE PRODUCT WHICH IS THE SUBJECT OF THE CLAIM UNDER WARRANTY. THE APPLICABILITY OF THE ABOVE MENTIONED LIMITATIONS AND EXCLUSIONS MAY VARY FROM STATE TO STATE AND/OR FROM

47

COUNTRY TO COUNTRY.

## SPECIFICATIONS OF HOWDY SENIOR

### Instructions for use:

Home/hospital monitoring of physiological parameters for adults over 16 years of age.

### Technical specifications of the electronic unit

Dimensions: 73x53x15 mm

Weight: 31gr.

Power supply: internal LiPo battery pack with a capacity of 600 mAh rechargeable with NFC technology (battery life 30 hours used with continuous data transmission through Bluetooth® connection).

NFC charger with standard Qi® (factor Q rx coil 25 – 125kHz)

Operating temperature range: +5°C - +45 °C

Operating relative humidity: 10% - 95% without condensation

Atmospheric operating pressure: 700-1060 hPa

Temperature range for storage/depot: -10 - +70°C

Relative humidity in storage/depot: 10% - 95% without condensation

Atmospheric pressure in storage/depot: 700-1060 hPa

Protection level of the electronic unit: IP22 (Protected against solid bodies of superior dimensions to 12mm, protected against falling water when tilted up to 15°)

Operating life without battery change: >500 measurements

Risk class: II b

Bluetooth version: 4.2  
NRF52832 module, transmission frequency band: 2,4 Ghz,  
Average consumption in transmission: 18.2 mA  
Characteristics of the modulation: Bluetooth v 4.2, Low Energy

#### ECG:

ECG channels: 1  
ECG recording: bipolar with two electrodes placed on the chest, and with a ground electrode below the left shoulder blade  
Frequency range: 0 - 512 Hz Resolution: 18 bits  
Heart rate calculation algorithm: Pan - Tompkins<sup>1</sup>.  
Automatic transmission via Bluetooth® 4.2 BLE and/or local storage of data

#### Breathing Signal

Acquisition through textile strain gauge integrated in the textile unit  
Frequency of sampling: 13Hz  
Resolution: 22 bit

#### Posture and Movement

Accelerometer: MEMS 3D  
Acceleration range:  $\pm 2g$   
Resolution: 15 bits

<sup>1</sup> Pan J, Tompkins W. A real-time QRS detection algorithm. *IEEE Trans Eng Biomed Eng.* 1985;32(3):230–236.

### **Electromagnetic Compatibility (EMC)**

#### **Manufacturer's declaration – Electromagnetic Emissions (Table 1, EN 60601-1-2)**

The device is suitable for use in a specific electromagnetic environment. The patient or family member using the device shall ensure that it is used in an electromagnetic environment as described below.

**Manufacturer's Guide and Declaration– electromagnetic emissions**

**The HOWDY SENIOR Device is intended to operate in the electromagnetic environment specified below. The customer or user should ensure that it is used in that environment**

| <b>Emission Test</b>     | <b>Compliance</b> | <b>Electromagnetic Environment</b>  |
|--------------------------|-------------------|---|
| Emissions RF<br>CISPR 11 | Group 1           | The medical device uses RF Energy only for its internal functioning. Therefore, the RF emissions are very low and avoids signal interferences with other nearby electronic devices.   |
| Emissions RF<br>CISPR 11 | Class B           | The medical device is suitable to be used in any environments, including domestic and all those environments that are directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |

**Manufacturer's declaration – Electromagnetic Immunity (Table 2, EN 60601-1-2)**

The device is suitable for use in a specific electromagnetic environment. The patient or family member using the device shall ensure that it is used in an electromagnetic environment as described below.

| <b>Immunity Test</b>                         | <b>IEC Test level 60601</b>                  | <b>Level of compliance</b>                   | <b>Electromagnetic Environment</b>   |
|--|--|--|--|
| Electrostatic discharge                      | ± 8 kV contact<br>± 2, 4, 8, 15 kV<br>in air | ± 8 kV contact<br>± 2, 4, 8, 15 kV<br>in air | Floor surfaces must be made of wood, concrete or ceramic tiles. If the floor surfaces are coated with synthetic materials, the relative humidity should be at least 30%. |
| Magnetic field at mains frequency (50/60 Hz) | 30A/m  | 30A/m  | Power frequency magnetic fields should be at levels characteristics of a typical location in a commercial or hospital environment.                                       |

**Manufacturer's declaration – Electromagnetic Immunity II (Table 4, EN 60601-1-2)**

The device is suitable for use in a specific electromagnetic environment. The patient or family member using the device shall ensure that it is used in an electromagnetic environment as described below.

| Immunity Test | IEC Test level 60601                          | Level of compliance                           | Electromagnetic Environment (distance between transmitter and device)  |
|---------------|---|---|--|
| Radiated RF   | 10 V/m<br>80 MHz – 2,7 GHz<br>80% AM to 1 kHz | 10 V/m<br>80 MHz – 2,7 GHz<br>80% AM to 1 kHz | $d = 0,35 \cdot \sqrt{P}$ from 80 MHz to 800 MHz<br>$d = 0,7 \cdot \sqrt{P}$ from 800 MHz to 1 GHz<br>$d = 2,33 \cdot \sqrt{P}$ from 1 GHz to 2,7 GHz<br>where P is the maximum output power rating in W according to the transmitter's manufacturer and represents the separation distance in metres. |

**Manufacturer's declaration - Recommended separation distances between portable and mobile RF communication equipment and the product (Table 206, EN 60601-1-2)**

The device is intended to operate in an electromagnetic environment where radio frequency radiated disturbances are under control. The patient or family member using the device can help prevent electromagnetic interference by ensuring a minimum distance between mobile and portable RF communication devices (transmitters) and the device as recommended below, in relation to the maximum output power of radio communication equipment.

| Potenza di uscita massima del Trasmittitore specificata (W) | Distanza di separazione alla frequenza del trasmettitore (m) |                   |                  |                  |
|---|--|-------------------|------------------|------------------|
|   | Da 15 kHz a 80MHz  | Da 80MHz a 800MHz | Da 800MHz a 1GHz | Da 1GHz a 2,5GHz |
| 0,01  | 0,12   | 0,04              | 0,07             | 0,23             |
| 0,1   | 0,37   | 0,11              | 0,22             | 0,74             |
| 1   | 1,17   | 0,35              | 0,70             | 2,33             |
| 10  | 3,69   | 1,11              | 2,21             | 7,38             |
| 100   | 11,67  | 3,50              | 7,00             | 23,33            |


### METHOD OF CALCULATING HEART RATE

The Pan and Tompkins algorithm is one of the most widely used methods for real-time recognition of the QRS complex from the ECG signal. The algorithm, based on slope, width and width of the QRS complex, includes a series of filters (low pass, high pass, shunt, quadrator and integrator) and methods (adaptive threshold and search procedure). The block diagram of the algorithm is shown in figure A.1.

In particular, the band-pass filter, obtained from the cascade of a low-pass filter and a high-pass filter, selects the frequency range in order to reduce the foreign components to the QRS, that is the P-wave, and the T-wave in particular; the shunt filter shows the rapid variation that characterizes the QRS complex; The quadrator filter makes the positive signal and emphasizes the components of the QRS complex and finally the moving average filter performs the filtering of the output signal from the previous filter, which could present multiple peaks at the QRS.

Figure a.1

*Block diagram of the Pan and Tompkins algorithm.*

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A large, empty rectangular box with a thin black border, occupying most of the page below the header text. It is intended for a receipt.